



Ace

SIEMENS



iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

Android and Google Play are trademarks of Google Inc. Apple App Store is a trademark of Apple Inc.

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks by Siemens AG is under license. Other trademarks and trade names are those of their respective owners.

The information in this brochure is subject to change without prior notice. © 2014 Siemens Hearing Instruments, Inc. All rights reserved.

Hearing instruments help many people hear better, but no hearing instrument can solve every hearing problem nor restore normal hearing.

Global Siemens Headquarters


Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Phone: +49 9131 84-0
Germany

usa.siemens.com/healthcare

usa.siemens.com/ace



touchControl App

Ace

Legal Manufacturer & Local Contact

Siemens Hearing Instruments, Inc.
P.O. Box 1397
Piscataway, NJ 08855-1397
usa.siemens.com/hearing

Global Business Unit

Siemens Audiologische Technik GmbH
Gebbertstrasse 125
DE-91058 Erlangen
Germany
Phone: +49 9131 308-0

Read the QR code with the QR code reader on your smart phone



usa.siemens.com/hearing

SH201461 7/14 1.0 22B771 1.5K SHI/15285-14

Enjoy the sound of life.

The new Ace™ hearing aid. Now more discreet than ever.



Life sounds brilliant.

The discreet duo for crystal-clear sound and convenient control. The new Ace with touchControl App.


As Siemens' tiniest Receiver-In-Canal (RIC) hearing aid, Ace™ offers you the utmost discretion and a superb listening experience. Plus, you can confidently and discreetly operate Ace with the new touchControl App: What looks like sending a text is really a hearing aid adjustment!

Discreet and powerful: RIC.



Unbelievably discreet, Ace's receiver sits directly in the ear canal – enabling a high level of amplification in an ultra-tiny housing.

At-a-glance: The new Ace.




A marvel of miniaturization: Although Ace is so tiny in size, it is packed with many sophisticated features. It is even resistant to moisture, sweat, dust and dirt according to IP67* standards.

- 1 touchControl option:**
Ace can be discreetly controlled via the touchControl App.
- 2 Directional microphones:**
Help understanding in demanding listening environments by focusing on sound from the front while reducing noise from other directions.
- 3 External receiver:**
With a choice of power levels, Ace is flexible to match your needs now and in the future.
- 4 Tinnitus therapy feature:**
A personalized sound signal helps to distract you from your tinnitus and may help to provide relief.
- 5 Push button:**
For manual program adjustments.
- 6 Battery compartment:**
Lets you easily turn your hearing aid on and off.



Ace's range of natural colors are designed to match a variety of hair and skin tones.



*Achieved IP67 rating per IEC 60529 standard. Water-resistant device can be completely submerged in water up to 3 feet for 30 minutes and dust will not interfere with the satisfactory operation of the device.

Ultimate discretion: With the new touchControl App.



Turn your smartphone into a discreet remote control for Ace and enjoy a high level of handling comfort right at your fingertips. The best part is: The download is app-solutely free!



touchControl App

With the touchControl App, you can easily and discreetly change Ace's hearing programs, adjust the volume and change the bass and treble.

You can even re-name the programs to make them more meaningful to you.

Easy to use and easy to get: App download.

In just a few simple steps, you can download the touchControl App on the Apple App Store or Google Play Store:

- ▶ Open your account on the Google Play Store or Apple App Store.
- ▶ Search for the touchControl App and download it to your Android or iOS smartphone.
- ▶ A brief on-screen guide will setup the touchControl App for your hearing aids.
- ▶ Enjoy full control right in your hands!



Enjoy life's true sound.

Ace is equipped with the latest version of BestSound™ Technology for the most natural hearing possible. It's a complete revolution in audiological engineering delivering a perfect balance of rich, natural sound quality and audibility. You may forget you're wearing hearing aids.

Ace is an optimum solution both for first-time and experienced wearers. Tiny and discreet in appearance, its enhanced BestSound Technology tailors the sound experience exactly the way you want to hear, right from the start.

To support the discretion of Ace, it can be controlled via an app for smart phones. The result: The latest BestSound Technology doesn't make you choose between discretion and functionality, you can now have both.

So wherever you are, whatever you do, Ace lets you enjoy the sound of life.



Life sounds brilliant.

